

BD- FACSCanto II

Lasers	Detection filter	Common detectable fluorochromes/ fluorescent proteins
Violet laser (405nm)	510/50	AmCyan; BV510; Qdot525
	450/50	BV421; Cascade Blue; Marina Blue; Pacific Blue; eFluor 450; DAPI; CFP; Alexa Fluor 405; Sytox Blue
Blue laser (488nm)	780/60	Pe-Cy7; Alexa Fluor750
	670LP	PE-Cy5; PerCP; PerCP-Cy5.5; PerCP-eFluor710
	585/42	PE; PI
	530/30	Alexa Fluor 488; FITC; GFP; CFSE; H2DCFDA; BB 515
Red laser (633nm)	780/60	APC-Cy7; APC-H7; APC-Alexa Fluor 750; APC-eFluor 780
	660/20	APC; Alexa Fluor 647; Alexa Fluor 633; Cy5; DDAO; eFluor 660; Sytox Red