


NMR ROOM 34



McGILL CHEMISTRY
CHARACTERIZATION

1. Room 34 contains: Varian Mercury 300 (unavailable as of May 21, 2020), **Bruker 400**, Varian 500 (unavailable as of May 21, 2020)
2. **Do not cut through room 34 to get to room 37**
3. Book time on FACES for the Bruker 400
4. Only two people are allowed in room 34 at a time:
 - One person waiting (sit at Varian 300)
 - One person using Bruker 400
 - Others must wait in the hall
5. Wear a mask
6. Gloves are optional
7. Wash your hands when you enter the room
8. Before doing NMR, you may use 95% ethanol and Kimwipes to wipe tube, spinner and depth gauge
9. Do not insert anything wet into the NMR
10. You may change the plastic wrap on the keyboard
11. After doing NMR, you must
 - Wipe with 70% ethanol: chair, table, handrail
 - Wipe with 95% ethanol: outside of spinner, depth gauge
 - Change plastic wrap on keyboard
12. Exit via room 35



Comments or suggestions? See <https://mcgill.ca/mc2/coronavirus>