


NMR ROOM 37



McGILL CHEMISTRY
CHARACTERIZATION

1. Room 37 contains: **Bruker 500**, QANUC 500 (unavailable as of May 21, 2020), **Varian 400**
2. **Do not cut through room 37 to get to room 34**
3. Book time on FACES for the Varian 400
4. Only two people are allowed in room 37 at a time (others wait in hall):
 - One person waiting (near the sink)
 - One person using Bruker 500
 - One person using Varian 400
5. Wear a mask
6. Gloves are optional
7. Wash your hands when you enter the room
8. Before doing NMR, you may use 95% ethanol and Kimwipes to wipe tube, spinner and depth gauge
9. Do not insert anything wet into the NMR
10. You may change the plastic wrap on the keyboard
11. After doing NMR, you must
 - Wipe with 70% ethanol: chair, table, handrail
 - Wipe with 95% ethanol: outside of spinner, depth gauge
 - Change plastic wrap on keyboard
12. Exit using same door you entered by



Comments or suggestions? See <https://mcgill.ca/mc2/coronavirus>